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INSTRUCTIONS FOR USE

Entnahmestellenblock Forano

Forano block of outlets

CE 0482

Deutsch 2

English 3

Entnahmestellenblock Forano

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Forano block of outlets

Dear customers,

thank you for purchasing this Greggersen product. If you have any questions or require information, please contact our sales and support team.

Please read the instructions for use thoroughly before using the medical device. Only use it once you have completely understood the application and function! All serious incidents must be reported to Greggersen Gasetechnik GmbH and to the health authority responsible for the place of residence.

Intended use

For connecting a maximum of three units to one outlet. Connection of units with plug-in nipples according to DIN 13260-2.

Operation

Before using the block of outlets, unplug all connected units from the plug-in couplings! It might be necessary to unplug other units plugged into the wall outlets. Plug the hose connection's angle plug into the outlet. Make sure to use the correct type of gas! For gas withdrawal, additional units can now be plugged into one of the outlets. After the first noticeable clicking, the plug-in nipple is retained in the parking position. The gas connection has not been established yet!

When the second snap-in position of the plug-in nipple is reached, the unit is in operating position: The gas connection has been established! The plug-in nipple can be removed from the terminal unit at any time by pressing on the actuator. The gas flow is interrupted automatically. In the process, hold the device / plug and press the actuator (single-hand operation). After use, unplug the hose connection's angle plug from the wall outlet.

Warnings

- The equipment and all connections **must** be free from oil, grease and lubricants since, otherwise, there is the risk of fire or explosion!
- Never use tools when connecting the unit and the accessories; all connections can be tightened or loosened by hand.

- Cleaning agents and disinfectants must not enter the unit!
- The total flow is limited by the tapping point which feeds the tapping point block.
- There is no warranty for correct function in case of improper transport or storage or during operating the block of outlets not in accordance to this instruction of use.

Maintenance / servicing

An external wiping disinfection is sufficient for all anodized aluminum parts. Clean plastic parts with a common disinfectant solution and follow manufacturers instructions. Completely dry the terminal unit prior to use!

Subject the terminal unit to a visual inspection prior to use (external damage). At least once a year, a function and leak test should be carried out. Servicing shall be carried out every six years, during which the cartridge has to be replaced (replacement of wearing parts). Faults and damage may be remedied only by authorised specialised staff. Use only original **GREGGERSEN** spare parts!

Spare parts and accessories

- 327463 Forano Cartridge DIN (other standards on request)
- Low-pressure-hose

Technical data

Operation pressure:	Compressed gas 400 ... 500 kPa, max. 1000 kPa, Vakuum \leq 40 kPa
Inlet:	NIST according to DIN EN ISO 18082
Outlet:	Plug socket according to DIN EN ISO 13260 part 2
Capacity:	up to 40l/min total flow
Weight:	2-fold 1020 g, 3-fold 1435 g
Temperature range:	-20 °C - 60 °C
Expected lifespan:	10 years
Guidelines:	DIN EN ISO 9170-1, DIN EN ISO 19054, DIN EN ISO 18082, DIN EN ISO 13260-2
Classification:	non-invasive; class IIa