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Instructions for use

Habicht

CE 0482

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
CE 0482

Intended use

For insufflation and extraction if central oxygen or medical air supply is not possible (emergencies, patient transfer, home-care, etc.).

Necessary accessories

- variant with Kolibri 2 or 3 l O₂ bottle
- variant with tube flow meter 3 l O₂ bottle
- pressure reducer
- low pressure hose

Operation

The instructions for use of the flow meters, suction units and pressure reducers apply.

Connect the pressure reducer to the O₂ bottle and attach it to the emergency supply unit at the NIST connection using a low-pressure hose line. The flow meter and the suction unit can now be put into operation. Remove pressure from the unit when replacing the cylinder, detach the hose from the pressure reducer and unscrew the connection nut of the pressure reducer from the cylinder valve by hand (no tools).

Warnings

- Prior to the use of the device, the user must be familiar with the product and the instructions for use and has to have understood them.
- The instructions for use of the flow meters, suction units and pressure reducers apply.
- The emergency supply unit Habicht with tube flow meter is not suitable for use with 2 liter bottles.
- Simultaneous use of both functions influences the performance and measurement accuracy.

- Never use tools when connecting the unit and the accessories; all connections can be tightened or loosened by hand. Always unscrew accessories when the pressure reducer is pressureless.
- The flowmeter will only display correctly if the measuring tube is vertical.
- Due to the oxidising effect of oxygen, special care must be taken when using oxygen-powered suction units.

Technical data


Type: Portable stainless steel frame, for suspension from an appliance rail or patient bed; with holder for gas cylinder of max. 3 liter (max. $\varnothing = 100\text{mm}$); Consisting of ejektor regulator and tube flowmeter or flowmeter Kolibri

Dimensions: approx. 500x330x150mm (HxWxD)

Medical Products Directive 93/42/ECC; Classification: IIa

Maintenance and inspection

The emergency supply unit must be serviced at least every 5 years. An inspection with a functional test and a leak test has to be performed at least every 2 years. In case of unusually heavy-duty use, shorter servicing and inspection intervals may be required. Servicing and repairs must only be carried out by authorised expert personnel. Only use genuine **Greggersen** spare parts.

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GREGGERSEN Gasetechnik GmbH / Bodestr. 27-29 / 21031 Hamburg / Germany