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
Kolibri 2

Dear customers,
thank you for purchasing this Greggersen product. If you have any questions or require information, please contact our sales and support team.

+49-(0)40 739 357-0, sales@greggersen.de

Preface

Please read the instructions for use thoroughly before using the medical device. Only use it once you have completely understood the application and function! The product Kolibri 2 from Greggersen Gasetechnik GmbH complies with the requirements of (EU) 2017/745 and is classified as a Class IIa active and non-invasive device. Basis UDI: 4031196KOLIBRIR7

Medical device in accordance with:

DIN EN ISO 15002

(EU) 2017/745

Intended use

Devices for administering medical gases.

Warnings

- The device may only be used for the intended purpose.
- Before each use, ensure the product is correctly aligned.
- Some of the products are made of brass with a lead content of over 0,1%. The lead is firmly bound in the material, so no exposure is to be expected.
- The product is not suitable for sterilization.
- No structural modifications may be made to the product.
- Nothing may be attached to the product that does not serve its direct use.
- Only connect or release NIST connections and accessories in the depressurised state. Never use tools for connecting the device or accessories. All connections can be released and tightened by hand.
- Oxygen is a highly fire-intensifying gas. Always keep all connections free of oil and grease! When handling, avoid the use of hand cream. Do not smoke or expose to open flame near oxygen equipment! Risk of fire!
- This product is intended exclusively for use in professional healthcare settings. Use in medical facilities is at the discretion of medical professionals.

- All serious incidents must be reported to Greggersen Gasetechnik GmbH and to the authority responsible for the place of residence.
- Ensure that the device does not generate any flow when in the intermediate position between two setting points

Variants

The devices are available in the following variants:

- Plug-in devices according to the DIN, BS, NF, SS, Medap and Carbamed standards
- Rail devices with NIST connection for DIN and Swiss standard device rails
- Different lengths available
- Single and double devices
- Devices for gas types O₂, AIR and CO₂
- other variants on request
- Devices with the following flow rate steps:

Flowlevel	l/min	l/min	l/min	l/min	l/min	l/min	l/min	l/min	l/min	l/min
Standard	0,5	1	1,5	2	3	4	6	8	12	15
Standard+	2	3	4	6	8	12	15	20	25	30
Children	0,1	0,3	0,5	0,6	0,8	1	2	3	4	5
Neonatal	0,1	0,15	0,2	0,25	0,3	0,4	0,5	0,6	0,8	1
o/8	8									

Marking

Consult the instructions
berücksichtigen

Do not use oil

Ordernumber

Temperaturerange

Medical device


Product contains lead

Manufacturer

(01) GTIN-Number
(11) Productiondate
(21) Serialnumber

Operation

Before each use, the product must be checked for proper functioning, mobility of the valves, as well as for external damage, obvious leaks and dirt. Non-functional, damaged, leaky or dirty devices must not be put into operation. Couple plug-in devices directly to a tapping point. Rail devices have to be attached to a device rail in accordance with DIN EN ISO 19054 and the low pressure hose has to be connected via the gas type coded NIST screw fitting. Before connecting, set the device to "0". Set the device with the hand wheel to the required value. Check whether the device generates an audible, tactile or (if a humidifier is used) visible gas flow. This ensures that a gas flow is actually generated and that the device is not in the parking position of the tapping point, for example. A connection to the patient can then be established either via a barbed connector at the 9/16-18 UNF connection or alternatively via a humidifier. The flow can be adjusted using the hand wheel to various values, which can be read through the viewing window.

After use, set the flow meter to "0". The device can then remain in the tapping point or be brought into the park position. If necessary, disconnect the humidifier or hose connection from the flow meter. The device can now be disconnected from the tapping point or released from the device rail. When using humidifiers, ensure that no liquid from the humidifier enters the device after it is removed from the tapping point.

Setting direction

Cleaning

It is recommended to clean the device by wiping disinfection after each use in accordance with the cleaning agent manufacturer's instructions. Commercially available disinfectants that are compatible with polycarbonate should be used. When cleaning, make sure that no liquid gets into the appliance. The device must be fully dried prior to reuse. Do not soak the device with disinfectant! The device may be damaged. The hygiene plan of the institution must be observed. The device is designed for multiple use with different patients.

Disposal

There are no special disposal requirements for this device. The device can be disposed of in household waste. If the product is contaminated, appropriate precautions must be taken.


Maintenance and inspection

Maintenance must be carried out every 5 years. An inspection including functional and leak testing must be performed at least every 2 years. Maintenance and repairs may only be carried out by Greggersen Gasetechnik or its authorized partners. The service life and maintenance interval begin on the date of commissioning, but no later than three years after the date of manufacture.

Device description

Pos	Description
1	Hand wheel
2	Viewing window with flow display
3	Gas type-specific plug connector
4	Rail clamp
5	Pressure inlet (NIST)
6	Flow outlet 9/16-18 UNF
7	Hose connection 9/16-18 UNF

Technical data

Gas supply:	Compressed air 400 to 500 kPa
Setting accuracy:	$\pm 20\%$; $\leq 1,5$ l/min $\pm 30\%$
Inlet:	Quick connector according to national standard (plug-in device) or NIST connection (rail-mounted device)
Outgoing connector:	9/16-18UNF External thread
Weight:	Single device Plug 340 g Double device plug 730 g Single unit rail 700 g Double device rail 1110 g
Dimensions:	See page 58
Storage conditions:	-20 to +60 °C; 15–93% relative humidity
Operating conditions:	+10 to +40 °C
Expected lifetime:	10 years

Spare parts and components necessary for operation

- 906470 Repair kit Kolibri 2 O₂
- 906471 Repair kit Kolibri 2 AIR
- 906485 Repair kit Kolibri 2 CO₂
- 900619 Hose connector 9/16-18UNF
- 904836 Humidifier 0,25 l
- Low pressure hoses

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